

TECHNICAL DATA



R7950/5L

Model R7950 Ultrasonic Couplant Gel, 250mL

- 8.45 fl oz (250ml) applicator top bottle
- Slow drying with good transducer lubrication
- Water soluble
- Includes applicator bottle
- Material Safety Datasheet (MSDS) available upon request

Model R7950/5L Ultrasonic Couplant Gel, 5L

- 5 Litre (1.32 gallon) barrel
- Slow drying with good transducer lubrication
- Water soluble
- Includes applicator bottle
- Material Safety Datasheet (MSDS) available upon request

Model R7950/12 Ultrasonic Couplant Gel, Pack of 12

- Slow drying with good transducer lubrication
- Water soluble
- Includes applicator bottle
- Material Safety Datasheet (MSDS) available upon request



R7950



R7950/12

| Model | Description |
|----------|-------------------------------------|
| R7950 | Ultrasonic Couplant Gel, 250mL |
| R7950/5L | Ultrasonic Couplant Gel, 5L |
| R7950/12 | Ultrasonic Couplant Gel, Pack of 12 |